


BRINGING YOU THE NEXT GENERATION OF IMPLANT EXTRACTION  
OSTEOTOMES



PIVOT  
IMPLANT EXTRACTION  
OSTEOTOME SYSTEM™



## 16 Different Modular Tips to Choose From:



Flat Straight blade short/medium/long/x-long



Curved blade short/medium/long



TM Peg Cutter



Straight Angled blade Short/medium/long



Angled Left and Right





The Patent Approved Pivot System™ is unique in that there is no comparable product designed specifically for removal of the knee component. The Pivot's modular osteotome tips take into account the complexity of the knee's different surface angles, flat areas and component sizes.

The unique design of the "V" shape, with three sharp sided blades and a pivot action point enables the surgeon to remove the knee components faster. The Pivot system reduces valuable OR time and minimizes bone loss.

The Pivot system has a wide variety of blade tips for the surgeon to choose from to work within the confines of the surgical area. Key benefits:

- Strong, thin blades
- Reduction in surgical time
- Preserves host bone
- Requires less impactful force to break the bond


## Pivot Design Features

### Ergonomically shaft design:

Giving the surgeon a perfect working distance

### Thin Blade:

This strong thin blade design is highly resistant to bending when using forceful impactions



### Pivot Mechanism:

Enables the surgeon to use the unique half circle mechanism to lever against all leading edges of the implant in a lateral left to right sweeping motion

### Three Sharp Sided-Blade:

Enables the surgeon to maneuver the blade in a lateral left to right motion with less force

### V-Shaped Blade:

Blade advances forward with less forceful impactions



## Pivot Osteotome Blades:

Reference No.	Description
39.1960.51.01	Right Angle, short
39.1960.51.02	Right Angle, long
39.1960.10.01	Straight, short
39.1960.10.02	Straight, medium
39.1960.10.03	Straight, long
39.1960.10.04	Straight, x-long
39.1960.20.01*	Straight Angled, short
39.1960.20.02	Straight Angled, medium
39.1960.20.03	Straight Angled, long
39.1960.30.01*	Curved, short
39.1960.30.02	Curved, medium
39.1960.30.03	Curved, long
39.1960.60.01*	Reverse Blade, short
39.1960.60.02	Reverse Blade, long
39.1960.52.01	Left Angle, short
39.1960.52.02	Left Angle, long
39.1960.70.01	Straight Wide, short
39.1960.70.02	Straight Wide, medium
39.1960.70.03	Straight Wide, long

## Pivot Osteotome Instruments:

Reference No.	Description
39.1960.01	Modular Handle
39.1960.02*	Tool Rod
39.1960.05	Retractor
39.1960.07	Slap Hammer
39.1960.100	Sterilization Case

\*Only Available in Canada

